


CERTIFICATES

We provide all the required documents such as manufacturing certificates, shipment documents and export documents for your tender requirements or projects quickly with our experience and expertise. CE, COO, ISO or other possible requirements.



Manufacturing, sales and after sales support for Dialysis equipment is available.



DIALYSIS CHAIR

A durable, high comfort, three electrical motor chair which can be used in Plasmapheresis treatment, has a special design and form. 360° rotatable wheels are with brakes and can be fixed to ground when required. All the metal parts are layered with electrostatic dye and anti-corrosive. The seats are semi-hard and covered with Finley leather which makes it easy to clean, antibacterial and water resistant.

Technical Specifications

Blood sampling chair is adjustable via electrical motor, you can adjust trendelenburg position with the control pad. There are two rails next to the chair and four knobs on the rail. Two armrests are synchronisedly functioning and moving with backseat position. You can dislocate and relocate armrests as per your requirement. Serum hanger is included in the package.

Dimensions

Width 560 mm - Length: 1850 mm - Height: 530 mm

Specifications Up and Down adjustable backseat and legs seats with remote control

All the surface is 5mm 28 Density gray rubber and covered with Vinlex leather.

Length as a patient chair 166 cm.

Width of back seat and feet seats 55 cm. Height from the ground 55 cm.

ExBxY – 55 cm x 166 cm x 55 cm



Designed with a wide range of movement options and functionality,

*** Colour options are available as per your requirement



REFERENCES

MEDICALPARK

ANADOLU^H

In Affiliation with
JOHNS HOPKINS MEDICINE

MEMORIAL



T.C. Sağlık Bakanlığı
Düzce Devlet Hastanesi



REFERANSANS

Azerbaijan Ministry of Health State Hospital, Dialysis Center technical service and spare part supply. Cameroon Ministry of Health Dialysis Centers consumables provider, Multiple Turnkey Dialysis center projects in Ghana.

Equitoreal Guinea Ministry of Health Dialysis Consumables supplier. Dialysis Center SetUp and Installment and long term supply partnership in Ukraine. Iraq State Hospital Water Treatment Reverse Osmosis System manufacturing and installment.

4 Kyrgyzstan Reverse Osmosis system and Solution manufacturing Plant



GALATA GLOBAL OUR DIALYSIS CENTERS



HD ISTANBUL DIALYSIS CENTER 35 Bads



UKRAYNA OMEGA DIALYSIS CENTER 15 Bads



GANA PSH DIALYSIS CENTER 22 Bads



THANKS

GALATA GLOBAL COMPANY

Ehlibeyt Mah. Tekstilciler Cad. Yapım İş Merkezi
No:35/14, Çankaya/Ankara

+90 312 247 50 25

info@galatamedical.com

www.galatamedical.com

